





exeter®

RETRIEVAL
SNARE

fingerspitzengefühl in fishing for foreign objects



- 1 Atraumatic Nitinol loop with gold-plated tungsten for highest visibility.
- 2 Pre-shaped MP Tip eases navigation.
- 3 Shapable introducer to ease navigation and gain performances.
- 4 AndraTec Supportive sleeve avoids kinking when manipulating.
- 5 Luer Connector.
- 6 Extra XXL torque device to ease manipulation
- 7 One piece kink-resistant wire enhances 1:1 torque for better manipulation.
- 8 Radiopaque platinum inner Marker Ring ensures nearest positioning of introducer.



The **Exeter Snare** and **Exeter Snare Micro** have been especially designed to meet the high requirements during extraction.

The **Exeter Retriever Technology** is your first choice safety device when it comes to successful retrieval.

Order codes

Product Code	⌀ (mm)	↓ (cm)	Supportive sleeve (cm/F)	Introducer sheath (F)
ES-05	5mm (3-5mm vessel)	125	110 / 4	4
ES-10	10mm (6-10mm vessel)	125	110 / 4	4
ES-15	15mm (11-15mm vessel)	125	110 / 4	4
ES-20	20mm (16-20mm vessel)	125	110 / 4	4
ES-25	25mm (21-25mm vessel)	125	110 / 5	5
ES-30	30mm (26-30mm vessel)	125	110 / 5	5
ES-35	35mm (31-35mm vessel)	125	110 / 5	5

Product Code	⌀ (mm)	↓ (cm)	Supportive sleeve (cm/F)	Introducer sheath (F)
ES-25/65	25mm (21-25mm vessel)	65	50 / 5	5
ES-25/150	25mm (21-25mm vessel)	150	135 / 5	5

Product Code	⌀ (mm)	↓ (cm)	Supportive sleeve (cm/F)	Introducer sheath (F)
ESM-2	2mm (1-2mm vessel)	175	150 / 3	3/4
ESM-4	4mm (3-4mm vessel)	175	150 / 3	3/4
ESM-7	7mm (5-7mm vessel)	175	150 / 3	3/4

Office & Warehouse

