



lokum

AMPLATZ TYPE
GUIDE WIRE

combining positioning and extra support in one design



- Design characteristics for positioning and extra support in one.
- High performance PTFE coated shaft providing low surface friction and enhanced trackability.
- Intended use during interventional peripheral and interventional biliary drainage procedures.
- Providing enhanced support during catheter and/or delivery system exchange.









Technical Information

| | |
|---------------------|------------------------------------|
| Guide Wire Material | Extra Tempered Stainless Steel |
| Coating | High performance PTFE |
| Diameter | 0.035" (0.89 mm) |
| Distal Tip | Straight / 3 mm J-Tip / 6 mm J-Tip |
| Tip Configuration | 5 cm Flexible Tip |
| Shaft Length | 180 cm / 260 cm / 360 cm |

Order Codes

A choice of various diameters and lengths.

|  | 3 mm  | 6 mm  |  |  |  |
|---|--|--|---|---|---|
| LOAM-180-ST | LOAM-180-3T | LOAM-180-6T | 0.035" (0.89 mm) | 180 cm | 5 cm |
| LOAM-260-ST | LOAM-260-3T | LOAM-260-6T | 0.035" (0.89 mm) | 260 cm | 5 cm |
| LOAM-360-ST | LOAM-360-3T | LOAM-360-6T | 0.035" (0.89 mm) | 360 cm | 5 cm |

Office & Warehouse

