




lokum

L-QUEST
GUIDE WIRE

balanced combination of extra strong support, high pushability
and variable tip lengths



- Especially designed for extra strong support high pushability and catheter / delivery system exchange.
- High performance PTFE coating providing low surface friction and enhanced trackability.
- Intended to use in complex diagnostic and interventional procedures.
- 3 tips flexibilities to ensure secure and safe positioning: 5 cm, 10 cm and 15 cm.



Technical Information

Lokum L-Quest Guide Wire







Guide Wire Material	Extra Tempered Stainless Steel
Coating	High performance PTFE
Diameter	0.035" (0.89 mm)
Distal Tip	Straight / 3 mm J-Tip / 6 mm J-Tip
Tip Configuration	5 cm / 10 cm / 15 cm Flexible Tip
Shaft Length	180 cm / 260 cm

Lokum L-Quest Mini Guide Wire







Guide Wire Material	Extra Tempered Stainless Steel
Coating	High Performance PTFE
Diameter	0.025" (0.635 mm)
Distal Tip	Straight / 3 mm J-Tip / 6 mm J-Tip
Tip Configuration	5 cm / 10 cm / 15 cm Flexible Tip
Shaft Length	180 cm / 260 cm

Order Codes

Lokum L-Quest Guide **0.035"**

	3 mm 	6 mm 			
LOLU-05-180-ST	LOLU-05-180-3T	LOLU-05-180-6T	0.035" (0.89 mm)	180 cm	5 cm
LOLU-10-180-ST	LOLU-10-180-3T	LOLU-10-180-6T	0.035" (0.89 mm)	180 cm	10 cm
LOLU-15-180-ST	LOLU-15-180-3T	LOLU-15-180-6T	0.035" (0.89 mm)	180 cm	15 cm
LOLU-05-260-ST	LOLU-05-260-3T	LOLU-05-260-6T	0.035" (0.89 mm)	260 cm	5 cm
LOLU-10-260-ST	LOLU-10-260-3T	LOLU-10-260-6T	0.035" (0.89 mm)	260 cm	10 cm
LOLU-15-260-ST	LOLU-15-260-3T	LOLU-15-260-6T	0.035" (0.89 mm)	260 cm	15 cm

Lokum L-Quest Mini Guide **0.025"**

	3 mm 	6 mm 			
LUMI-05-180-ST	LUMI-05-180-3T	LUMI-05-180-6T	0.025" (0.635 mm)	180 cm	5 cm
LUMI-10-180-ST	LUMI-10-180-3T	LUMI-10-180-6T	0.025" (0.635 mm)	180 cm	10 cm
LUMI-15-180-ST	LUMI-15-180-3T	LUMI-15-180-6T	0.025" (0.635 mm)	180 cm	15 cm
LUMI-05-260-ST	LUMI-05-260-3T	LUMI-05-260-6T	0.025" (0.635 mm)	260 cm	5 cm
LUMI-10-260-ST	LUMI-10-260-3T	LUMI-10-260-6T	0.025" (0.635 mm)	260 cm	10 cm
LUMI-15-260-ST	LUMI-15-260-3T	LUMI-15-260-6T	0.025" (0.635 mm)	260 cm	15 cm

Office & Warehouse

